Yes, you can.®

EN

Power Wheelchair Wiring Guide

NOTICE

This manual contains recommended guidelines for the proper routing and securing of wires, cables, and tubing when removing/installing or replacing electronic components on Invacare power wheelchairs or conducting periodic inspection of the wheelchair.

The following are guidelines **ONLY**. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories.

Decisions regarding routing and securing of wires and cables should **ONLY** be made by a qualified technician.

DANGER

Risk of Death, Serious Injury, or Damage

Incorrect repair and/or servicing of this wheelchair performed by users/caregivers or unqualified technicians can result in death, serious injury, or damage.

- Users/Caregivers — **DO NOT** attempt to repair and/or service this wheelchair.
- Repair and/or service of this wheelchair **MUST** be performed by a qualified technician. Contact a dealer or Invacare technician.

WARNING

Risk of Death, Injury or Damage

Improper use of this product may cause injury or damage

If you are unable to understand the warnings, cautions or instructions, contact a health care professional or dealer before attempting to use this equipment.

DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manuals or instruction sheets supplied with this product or optional equipment.

Risk of Injury or Damage

Incorrect set up of this wheelchair performed by users/caregivers or unqualified technicians can result in injury or damage.

User/Caregivers- **DO NOT** attempt to set up this wheelchair.

Initial set up of this wheelchair **MUST** be performed by a qualified technician.

Service Manual

DEALER: Keep this manual. The procedures in this manual **MUST** be performed by a qualified technician.

TABLE OF CONTENTS

SECTION 1—SAFETY 3

- Symbols.....3
- Warnings.....3
- General Guidelines3
- Service Information.....4
- G-Trac.....4
- Sip n’ Puff (Pneumatic) Controls5
- Grounding Instructions5
- Batteries.....5
- Cleaning6

SECTION 2—WIRING 7

- Wire Routing Guidelines7
- Common Causes of Wiring Damage.....7
- Potential Results of Wiring Damage7
- Pinch Points.....7
- Removing/Installing or Replacing Electronic Components9
- Bundling the Wiring..... 10
- Recommended Wire Routing and Securing Points 11
- Wheelchair Base 11
- Adjustable ASBA and Van Seat Arms 12
- Conventional Back..... 13
- Contoura Back 14
- Powered Seating 15

I Safety

I.1 Symbols

Warnings

Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

DANGER

Danger indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Gives useful tips, recommendations and information for efficient, trouble-free use.

I.1 General Guidelines

The Safety section contains important information for the safe operation and use of this product.

WARNING

Risk of Death, Injury or Damage

Improper use of this product may cause injury or damage

If you are unable to understand the warnings, cautions or instructions, contact a health care professional or dealer before attempting to use this equipment.

DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manuals or instruction sheets supplied with this product or optional equipment.

Risk of injury or damage

Incorrect repair and/or servicing of this wheelchair performed by users/caregivers or unqualified technicians can result in serious injury or damage.

Users/Caregivers — DO NOT attempt to repair and/or service this wheelchair.

Repair and/or service of this wheelchair MUST be performed by a qualified technician. Contact a dealer or Invacare technician.

ACCESSORIES WARNINGS

Risk of injury or damage

Use of non-Invacare accessories may result in injury or damage.

Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories. Accessories designed by other manufacturers have not been tested by Invacare and are not recommended for use with Invacare products.

DO NOT use non-Invacare accessories.

To obtain Invacare accessories, contact Invacare by phone or at www.invacare.com

DO NOT connect any medical devices such as ventilators, life support machines, etc., directly to the batteries used to power the wheelchair. This could cause unexpected failure of the device and the wheelchair.

Wiring must never be modified or altered in any way to provide power tapes, connections or features not originally part of the wheelchair such as chargers for cell phones, music players, game systems, lap top computers, lighting or other electrical devices.

REPLACEMENT PARTS WARNINGS

Risk of injury or damage

Use of incorrect or improper replacement (service) parts may cause injury or damage.

Replacement parts **MUST** match original Invacare parts.

Always provide the wheelchair serial number to assist in ordering the correct service parts.

NOTICE

THE INFORMATION CONTAINED IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE.

I.2 Service Information

DANGER

Risk of Death, Serious Injury, or Damage

Continued use of the wheelchair that is not set to the correct specifications may cause erratic behavior of the wheelchair resulting in death, serious injury, or damage.

Performance adjustments should only be made by professionals of the healthcare field or persons fully conversant with this process and the driver's capabilities.

After the wheelchair has been set up/adjusted, check to make sure that the wheelchair performs to the specifications entered during the set up procedure. If the wheelchair does not perform to specifications, turn the wheelchair Off immediately and reenter set up specifications. Contact Invacare, if wheelchair still does not perform to correct specifications.

WARNING

Risk of injury or damage

Set-up of the driver control is to be performed only by a qualified technician. The final adjustments of the driver control may affect other activities of the wheelchair. Damage to the equipment could occur if improperly set-up or adjusted.

Before adjusting, repairing or servicing the wheelchair, **ALWAYS** turn the wheelchair power Off, otherwise, injury or damage may occur.

Except for programming, **DO NOT** service or adjust the wheelchair while occupied, unless otherwise noted.

Wheelchairs should be examined during maintenance for signs of corrosion (water exposure, incontinence, etc.). Electrical components damaged by corrosion should be replaced immediately.

Wheelchairs that are used by incontinent users and/or are frequently exposed to water may require replacement of electrical components more frequently.

DO NOT overtighten hardware attaching to the frame. This could cause damage to the frame tubing.

After **ANY** adjustments, repair or service and before use, make sure all attaching hardware is tightened securely - otherwise injury or damage may occur.

G-Trac

WARNING

Risk of injury or damage

Improperly connecting the motor leads to the controller may cause injury or damage.

WHEELCHAIRS WITH G-TRAC: Crossing the motor leads (for example: connecting the left motor lead into the right motor connector on the controller) may result in unintended movement.

DO NOT cross the motor leads when connecting the motors to the controller - otherwise injury or damage may occur.

Sip n' Puff (Pneumatic) Controls

WARNING

Risk of injury or damage

Improper mounting or maintenance of the Sip n' Puff control including the mouthpiece and breath tube may cause injury or damage.

Water inside the Sip n' Puff interface module may cause damage to the unit.

Excessive saliva residue in the mouthpiece/straw can reduce performance.

Blockages, a clogged saliva trap or air leaks in the system may cause Sip n' Puff not to function properly.

Ensure moving parts of the wheelchair, including the operation of powered seating, **DO NOT** pinch or damage the Sip n' Puff tubing.

Saliva trap **MUST** be installed to reduce risk of water or saliva entering the Sip n' Puff interface module.

Occasionally flush the mouthpiece to remove saliva residue.

The mouthpiece/straw **MUST** be completely dry before installation.

If Sip n' Puff does not to function properly, inspect system for blockages, clogged saliva trap or air leaks.

As necessary, replace mouthpiece, breath tube and saliva trap.

Grounding Instructions

WARNING

Risk of injury or damage

Improper grounding or use of three prong adapters may cause injury or damage.

DO NOT, under any circumstances, cut or remove the round grounding prong from any plug used with or for Invacare products. Some devices are equipped with three-prong (grounding) plugs for protection against possible shock hazards and fire. Where a two prong wall receptacle is encountered, it is the personal responsibility and obligation of the customer to contact a qualified electrician and have the two prong receptacle replaced with a properly grounded three prong wall receptacle in accordance with the National Electrical Code. If you must use an extension cord, use **ONLY** a three-wire extension cord having the same or higher electrical rating as the device being connected. In addition, Invacare has placed **RED/ORANGE WARNING TAGS** on some equipment. **DO NOT** remove these tags.

Three prong to two prong adapters should not be used. Use of three prong adapters can result in improper grounding and present a shock hazard to the user.

Batteries

WARNING

Risk of injury or damage

Improper use or care of batteries may cause injury or damage.

The use of rubber gloves is recommended when working with batteries.

After **ANY** adjustments, repair or service and before use, make sure all attaching hardware is tightened securely - otherwise injury or damage may result.

NEVER allow any of your tools and/or battery cable(s) to contact **BOTH** battery post(s) at the same time. An electrical short may occur and serious personal injury or damage may occur.

The **POSITIVE (+) RED** battery cable **MUST** connect to the **POSITIVE (+)** battery terminal(s)/post(s), otherwise serious damage will occur to the electrical system.

Install protective caps on **POSITIVE (+)** and **NEGATIVE (-)** battery terminals.

DO NOT remove fuse or mounting hardware from **POSITIVE (+) RED** battery cable mounting screw.

Cleaning

WARNING

Risk of damage

Improper cleaning may cause damage.

DO NOT use spray cleaning chemicals, polishes, waxes, oils or other similar products on wiring, connectors or electrical components of the wheelchair. These chemicals can degrade the insulation properties or accelerate wear.

Cleaning should only be performed with a water damped towel.

NEVER use water forced from a garden hose or power washer to clean the wheelchair.

Updated versions of this manual are available on www.invacare.com

2 Wiring

The following are guidelines ONLY. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories.

Decisions regarding routing and securing of wires and cables should only be made by a qualified technician.

If technical questions or problems arise, contact Invacare Technical Service.

In this manual the term “wiring” describes any cable, lead or wire used to connect a electronic/powering device to the wheelchair batteries or charger.

2.1 Wire Routing Guidelines

It is important to follow these basic wiring guidelines when installing, adding, repositioning or replacing wiring of electronic components such as the controller, motors, joystick, display, switches or accessory/seating control box(es). This manual points out common procedures and techniques involved in the safe installation, repositioning or maintenance of power wheelchair wiring.

Use this information only as a “basic” guide. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories.

Installation, repositioning or maintenance of power wheelchair wiring requires close attention to ensure wiring does not interfere with the safe operation of the wheelchair and to prevent damage to the wiring, the wheelchair and/or injury to the user.

Service personnel must be aware that the location and wire routing and location of mounting hardware can effect the following:

- User safety.
- Performance of the wheelchair or accessories.
- Location of driver controls such as joysticks, switches and displays.
- The range of motion of accessories such as powered seating and front riggings.

Great care must be used to ensure all wiring is properly routed and secured to the wheelchair to minimize risk of pinched or damaged wiring and to prevent loose wiring from snagging on objects around the wheelchair.

Before attempting to reposition, remove or replace any existing wiring, note the position and orientation of the following:

- The wire routing along the wheelchair frame.
- Location of tie-wraps and any other hardware securing the wiring to the wheelchair.
- location of wiring connections.

Wiring should only be repositioned if necessary to make an adjustment to fit the wheelchair to the user or if needed to add a new option - Otherwise wiring must be reinstalled, routed and secured in the original path along the wheelchair frame.

2.2 Common Causes of Wiring Damage

Improperly routed wiring may become damaged under the following conditions:

- Wires left loose or improperly secured.
- Wires secured too tightly to the wheelchair.
- Wires routed through pinch points.
- Wiring improperly secured to moving parts.

Ensure wiring has been properly secured to the wheelchair and has not been left loose or hanging from the wheelchair. Otherwise wiring could be snagged or entangled during normal operation of wheelchair. Cables shall be inspected periodically to ensure that they are routed and secured properly. Periodic inspection is recommended as it may reveal loose and/or damaged cables. Re-secure all loose cables and replace all damaged cables by following the recommendations outlined in this manual.

Moving parts of a wheelchair such as powered seating systems, wheelchair suspension and side frame/swing arm may cause a pinch point where wiring could become cut or crushed resulting in damage to the wire, wheelchair or user. Avoid routing wires through or near pinch points.

Wiring may also be damaged if improperly routed through areas of repetitive motion such as a quadlink or elevating seat. Ensure wiring in these areas allow for full range of motion without leaving excess wiring that could become snagged or entangled while driving.

Potential Results of Wiring Damage

Failure to follow the guidelines to properly install accessories, or to route wiring may lead to:


- Poor performance of wheelchair or accessories.
- Electrical fires.
- Intermittent electrical arcing or shorting.
- Personal injury or death of user or others.
- Potential for smoking, melting or damage to wiring.

Pinch Points


WARNING

A pinch point exists between the headrest knob/bracket and recline actuator. When routing the wiring, ensure that all wires are clear of pinch point area. Otherwise, serious injury or damage may occur. See Detail “A”.


CORRECT WIRE ROUTING

i Refer to Bundling the Wiring on page 10 for correct routing of excess wires.



For views of wiring for Conventional Back, refer to FIGURE 2.5 on page 13. For views of wiring for Contoura Back, refer to FIGURE 2.6 on page 14.

INCORRECT WIRE ROUTING

Wires should NOT be routed in this manner as they may get pinched or severed during the recline function.

i Examples of common pinch points and areas to avoid when securing wiring to the wheelchair.
DO NOT wrap wiring around moving parts or secure in a location where wires maybe crushed or damaged.

i DO NOT secure wiring to rear suspension springs or insert wiring through the rear suspension springs.


FIGURE 2.1 Common Causes of Wiring Damage

2.3 Removing/Installing or Replacing Electronic Components

The following are recommended guidelines for the routing and securing of wiring while removing/installing or replacing electronic components or wiring on Invacare power wheelchairs.

The following are guidelines **ONLY**. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories.

1. Before installing any electrical component or mounting hardware, including tie-wraps, cable clamps or mounting brackets, perform the following to ensure proper installation of the electrical component(s):
 - Note position and orientation of tie-wraps securing the wiring to the wheelchair.
 - Note the wiring path along the wheelchair frame.
 - If removing tie-wraps securing several wires together or wires that have been folded together as a bundle, note the wires included in bundle and how the wires have been folded together.
 - When removing an electrical component such as the controller, accessory control box, switch, display or joystick note position and orientation of mounting hardware.
2. When installing new/existing electronic components ensure all components and wiring are replaced in the original position and orientation unless it is necessary to reposition the components to allow for installing a new component or to adjust the wheelchair to the user.
3. **After adjustment or service perform the following checklist:**
 - Inspect all fasteners.**
 - Ensure tie-wraps are tight enough to prevent movement of wires under normal use without deforming the wiring insulation.**
 - Ensure wiring connections are secure.**
 - Inspect electrical components for signs of corrosion. Replace if corroded or damaged.**
 - Verify charger function.**
 - Inspect wiring attached to moving parts, such as retractable arms or seating systems for signs of chaffing, abrasion, stretching or for sections of wire that appear visually different than the other wires in the same area.**
 - With the wheelchair unoccupied, test all system functions, multiple times, through full range of motion to verify proper operation and to verify wiring harnesses do not obstruct the path of the system. If they do, perform one of the following:**
 - **Wires were damaged during inspection - Replace damaged wires ensuring these are properly located and secured where they will not become damaged.**
 - **Wires were not damaged during inspection - Cut tie-wraps and relocate wires to a location where they will not become damaged.**
4. If technical questions or problems arise, contact Invacare Technical Service.

2.4 Bundling the Wiring

i In this manual, the terms bundle and bundling describe a method of folding and binding excess wiring together to reduce risk of damage from loose or improperly secured wiring.

In some situations, it may be necessary to bundle excess wiring before securing the wiring to the wheelchair frame. Refer to the following guidelines for bundling wiring:

1. Starting at one end of the cable, secure the wire to the wheelchair.
2. Alternating ends of the cable, continue securing the wire to the frame.
3. If necessary, gather excess wiring and fold the wire together into as compact a bundle as possible.

i Before binding bundle together, ensure there is an appropriate position on the wheelchair to secure the bundle. If necessary, release the bundle and perform STEP 3 to change the bundle location.

4. Using one to two tie-wraps, bind the bundle together.
5. Using one to two tie-wraps, secure the bundle to the wheelchair.

i Examples of excess wiring gathered into bundles before being secured to the wheelchair. Images are examples of wire bundling ONLY. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories. ALL WIRING, CONNECTOR STACK, AND CONTROL BOXES **MUST** BE PROPERLY SECURED TO THE WHEELCHAIR BEFORE USE.



FIGURE 2.2 Bundling the Wiring



2.5 Recommended Wire Routing and Securing Points

The following figures contain examples of common wire routing and securement points, however, it must be noted that these examples are presented as guidelines ONLY.


- i** The following are guidelines ONLY. Location of electronic components and wiring may vary from wheelchair to wheelchair according to wheelchair model, manufacturing date, options or accessories.
Decisions regarding routing and securing of wires and cables should only be made by a qualified technician.
If technical questions or problems arise, contact Invacare Technical Service.

Wheelchair Base

- i** Examples of typical tie-wrap locations built into the rear suspension and frame of TDX SP, SR, SI, SC, and Spree Wheelchairs.

- i** Example of typical tie-wrap locations on the rear suspension pivot tube of TDX™ SP, SR, SI, SC, and Spree Wheelchairs.


- i** Example of typical tie-wrap locations on the rear battery bracket of the FDX™ .


FIGURE 2.3 Recommended Wire Routing and Securing Points - Wheelchair Base

Adjustable ASBA and Van Seat Arms

i Typical tie-wrap locations on joystick tube without a quadlink joystick mount.

i Examples of typical tie-wrap locations on joystick tube with a quadlink joystick mount.


i The use of a cable bracket with a quadlink joystick mount is strongly recommended to reduce stress on the joystick cable.

FIGURE 2.4 Recommended Wire Routing and Securing Points - Adjustable ASBA and Van Seat Arms

Conventional Back

i Example control boxes and wiring secured to conventional back.

i Examples of wiring locations on joystick/display mounting tubes.



FIGURE 2.5 Recommended Wire Routing and Securing Points - Conventional Back

Contoura Back

Example control boxes and wiring secured to Contoura back.

Examples of wiring locations on joystick tubes/display/mounting tubes.



FIGURE 2.6 Recommended Wire Routing and Securing Points - Contoura Back

Powered Seating

i Example powered seating wiring secured to wheelchair frame.


FIGURE 2.7 Recommended Wire Routing and Securing Points - Powered Seating

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© 2014 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. Phillips is a registered trademark of Phillips Screw Company. Teflon is a registered trademark of DuPont.

I167603-D-05

2014-12-9

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Invacare Corporation www.invacare.com